



FlowWeave[®] plus
Blood-Tight Woven Polyester Grafts

FlowWeave[®] plus

Blood-Tight Woven Polyester Grafts

Thanks to a special weaving technique, FlowWeave plus vascular grafts have properties which minimize porosity giving excellent haemostatic results.


The use of texturized yarn provides for superior tissue incorporation and healing characteristics.

JOTEC's modern weaving technique guarantees high burst strength and low dilatation.


Excellent handling and ease of suturing of FlowWeave plus vascular grafts facilitates operating techniques, and helps reduce intervention times.



FlowWeave plus tetra-branch thoracic with integrated side branch for antegrade perfusion minimizes circulatory arrest time and reduces risk of neurological deficits



Reconstruction of aortic arch and supra-aortic vessels




Reconstruction of thoracoabdominal aorta along with efferent major arteries

FlowWeave plus tetrabranch abdominal with four side branches enables sequential reconnection of arteries and restoration of blood flow

FlowWeave plus vascular grafts with and without side-branches are indicated for reconstruction of thoracic and abdominal aorta especially in case of stenoses, aneurysms and dissections.

Blood-tight without impregnation and preclotting – FlowWeave plus



Concentric crimping together with guide lines enables safe and exact positioning

FlowWeave[®] plus

Blood-Tight Woven Polyester Grafts

Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49 (0) 7471 / 922-0

FlowWeave[®] plus

Length (cm)	ø (mm)	Catalogue-No.	Length (cm)	ø (mm)	Catalogue-No.
15	24	41ST1524	30	24	41ST3024
15	26	41ST1526	30	26	41ST3026
15	28	41ST1528	30	28	41ST3028
15	30	41ST1530	30	30	41ST3030
15	32	41ST1532	30	32	41ST3032
15	34	41ST1534	30	34	41ST3034
15	36	41ST1536	30	36	41ST3036
15	38	41ST1538	30	38	41ST3038

FlowWeave[®] plus tetrabranch thoracic

Usable Length Main Tube 50 cm, Side Branches 15 cm

ø Main Tube (mm)	ø Branch 1 (mm)	ø Branch 2 (mm)	ø Branch 3 (mm)	ø Branch 4 (mm)	Catalogue-No.
20	10	8	8	10	41BR2010080810T
22	10	8	8	10	41BR2210080810T
24	10	8	8	10	41BR2410080810T
26	10	8	8	10	41BR2610080810T
28	10	8	8	10	41BR2810080810T
30	10	8	8	10	41BR3010080810T
32	10	8	8	10	41BR3210080810T
34	12	10	10	10	41BR3412101010T
36	12	10	10	10	41BR3612101010T
38	12	10	10	10	41BR3812101010T

FlowWeave[®] plus tetrabranch abdominal

Usable Length Main Tube 50 cm, Side Branches 15 cm

ø Main Tube (mm)	ø Branch 1 (mm)	ø Branch 2 (mm)	ø Branch 3 (mm)	ø Branch 4 (mm)	Catalogue-No.
20	10	10	8	8	41BR2010100808A
22	10	10	8	8	41BR2210100808A
24	10	10	8	8	41BR2410100808A
26	10	10	8	8	41BR2610100808A
28	10	10	8	8	41BR2810100808A
30	10	10	8	8	41BR3010100808A
32	10	10	8	8	41BR3210100808A

JOTEC[®]

SOLUTIONS FOR VASCULAR DISEASE

JOTEC AG

Luzernerstr. 91
CH-5630 Muri
T +41 (0) 56 670 91 - 00
F +41 (0) 56 670 91 - 02
www.jotec.com

JOTEC GmbH

Lotzenerstr. 23
D-72379 Hechingen
T +49 (0) 74 71 / 922 - 0
F +49 (0) 74 71 / 922 - 100

JOTEC s.r.l.

Vicolo Valtellina n.15
I-20092 Cinisello Balsamo (Mi)
T +39 02 30310061
F +39 02 66014455

JOTEC CARDIOVASCULAR SL

Av. Castilla 2, P. E. San Fernando
Edificio Italia – planta baja
E-28830 San Fernando de Henares, Madrid
T +34 (0) 91 - 648 70 20
F +34 (0) 91 - 656 17 73

